


Vivaspin® 4 and Vivaspin® 6

Concentration & Purification of Clinical Biological Samples



Vivaspin® ultrafiltration concentrators are disposable devices for biological samples. Vivaspin® 4 is suitable for sample volumes up to 4 ml while Vivaspin® 6 can handle samples as much as 6 ml. The patented vertical membrane design and thin channel filtration chamber (US 5,647,990) minimizes membrane fouling and provides high speed concentrations. Vivaspin® devices are for single use and Vivaspin® 4 devices are intended for in vitro diagnostic (IVD) applications. Some Vivaspin® 6 products are also intended for IVD applications.

Vivaspin® Applications

- Concentration of urine prior to electrophoresis for diagnosis of multiple myeloma and amyloidosis.
- Concentration of spinal fluid prior to electrophoresis for diagnosis of meningitis and multiple sclerosis.
- Concentration of bacterial antigens in urine (*Legionella*, *Pneumonia*, *Streptococcus B*) prior to immunoassay

Features	Benefits
Vertical membrane design	Fast filtration rates.
Integrated dead stop	Reduces risk of concentration to dryness.
Low binding materials	High sample recovery.
Wide range of MW cutoffs	Choose MW cutoff for your specific application.
Large Vivaspin® family of devices	Can process samples from 100 µl to 20 ml with similar methodology.

Technical Specifications

Concentrator Capacity	Vivaspin® 4	Vivaspin® 6
Fixed angle rotor	4 ml	6 ml
Swinging bucket rotor	4 ml	6 ml

Dimensions

Total Length	122 mm	122 mm
Width	17 mm	17 mm
Active membrane area	2.0 cm ²	2.5 cm ²
Hold-up volume	< 10 µl	< 10 µl
(membrane & support)		
Dead stop volume**	20 µl	30 µl
**Dead stop volume may vary depending on sample, temperature and rotor		

Materials of Construction

Body	Polycarb.	Polycarb.
Filtrate vessel	Polyprop.	Polycarb.
Concentrator cap	Polycarb.	Polyprop.
Membrane *	PES	PES

Polycarb. = Polycarbonate, Polyprop. = Polypropylene
* See other side for membrane descriptions

Equipment Required

Centrifuge

Rotor cavity	To fit 15 ml (17 mm) conical bottom tubes	To fit 15 ml (17 mm) conical bottom tubes
Fixed Rotor - Minimum Angle	25°	25°
Fixed Rotor - Max. Speed	8,000 g	8,000 g **
**lower g force for these membranes: ≥ 100,000 PES	N/A	6,000 g
Swinging Bucket -Max. Speed	4,000 g	4,000 g

Pipettes for Sample Delivery & Recovery

Fixed or variable volume may be used.
For maximum recovery, a thin gel loader type is recommended.

Membrane Selection Guide

Polyethersulfone (PES)

A good general purpose membrane for most solutions. Has low fouling characteristics with very good flow rates. Tolerates a broad pH range (from 1 to 9).

Cellulose Triacetate (CTA)

Good hydrophilic properties with very low non-specific binding. Cast without a membrane support that could bind filtered solutes. Use when recovery of filtrate is most important (free drug or hormone testing).

Hydrosart® (HY)

Similar to regenerated cellulose but with better performance characteristics and low protein binding. Good choice for concentration and desalting of proteins that may bind to other membranes such as immunoglobulin fractions.

Performance Characteristics

Time (in min.) to concentrate 30x at 20° C with rotor type shown. Also showing sample % recovery.

	Vivaspin® 4 4 ml@5,000 g Fixed angle		Vivaspin® 6 6 ml@7,500 g Fixed angle		Vivaspin® 6 6 ml@3,000 g Swinging Bucket	
	Time	Rec.	Time	Rec.	Time	Rec.
Cytochrome c-0.25 mg/ml (12,400 MW) 5,000 MWCO PES	N/A		90	97%	N/A	
BSA-1.0 mg/ml (66,000 MW)						
5,000 MWCO PES	N/A		12	98%	20	98%
10,000 MWCO PES	10	96%	10	98%	13	98%
30,000 MWCO PES	N/A		9	97%	12	98%
IgG-0.25 mg/ml (160,000 MW)						
30,000 MWCO PES	N/A		15	95%	18	96%
50,000 MWCO PES	N/A		14	95%	17	96%
100,000 MWCO PES	N/A		12	91%	15	91%

Ordering Information

Molecular Weight Cutoff (MWCO)	Pack Size	Vivaspin® 4 PES membrane	Vivaspin® 6 PES membrane
3,000 MWCO	25		VS0691
3,000 MWCO	100		VS0692
5,000 MWCO	25		VS0611
5,000 MWCO	100		VS0612
10,000 MWCO	25	VS0403	VS0601IVD
10,000 MWCO	100	VS0404	VS0602IVD
10,000 MWCO	500	VS0405	
30,000 MWCO	25		VS0621
30,000 MWCO	100		VS0622
50,000 MWCO	25		VS0631
50,000 MWCO	100		VS0632
100,000 MWCO	25		VS0641
100,000 MWCO	100		VS0642
300,000 MWCO	25		VS0651
300,000 MWCO	100		VS0652
1,000,000 MWCO	25		VS0661
1,000,000 MWCO	100		VS0662
0.2 µm	25		VS0671
0.2 µm	100		VS0672

Manufactured by: Sartorius Stedim Lab Ltd. • Sperry Way, Stonehouse Park
GL10 3UT Stonehouse, Gloucestershire, UK • www.sartorius.com

North American Distributors:



521 Great Road • Littleton, MA 01460 • www.vivaproducts.com • Phone: 800.456.4633 • Fax: 877.289.8482

Visit our website www.vivaproducts.com for pricing and online ordering!