



Application Guide

December 05, 2017

Sartorius Ultrafiltration Products in the Preparation of Biological Nanoparticles and Medical Nanocarriers – a Short Review

Hannes Landmann^{1*}, Kristin Menzel²

¹ Sartorius Lab Instruments GmbH & Co. KG, Otto-Brenner-Straße 20, 37079 Göttingen, Germany

² Sartorius Stedim Biotech GmbH, August-Spindler-Straße 11, 37079 Goettingen, Germany

* Correspondence

E-Mail: hannes.landmann@sartorius.com

Abstract

Keywords or phrases:

Nanoparticles,
Nanocarriers,
Liposomes, Vesicles,
Micelles

This short review highlights ultrafiltration of various biological nanoparticles and medical nanocarriers. Here ultrafiltration is used to purify, concentrate and separate nanoparticles from substrates. The nanomaterials discussed comprise metals, polymers, lipids (in vesicles) and proteins. You will find a guidance for the selection of an ideal performing ultrafiltration device with the optimum molecular weight cut-off (MWCO) for typical concentration applications.

Introduction

Paul Ehrlich was inspired by the idea of the "magic bullet"* when he for the first time described in theory toxic drugs assembled to so-called "Nanocarriers" in 1908.¹ Today, Nanocarriers have found multiple applications in modern medicine and biotechnology. A key application for these special nanomaterials is a targeted delivery of drugs where they act as transport modules (i. e. as nanoparticles, vesicles, or micelles) for the active ingredient.^{2,3,4,5} This is assumed to be more effective and less toxic to the (human) organism compared to traditionally administered drug substances.⁶ Besides drug delivery, various further fields using Nanocarriers evolved during the last decades; e. g. magnetic resonance imaging or stem cell gene therapy with metal-based nanoparticles,^{7,8} or optical imaging with quantum dots.⁹

Nanocarriers can be categorized by their starting material (i. e. metal-, lipid-, polymer-, and protein-based) and by their formation after preparation (i. e. vesicles, particles and micelles). In general, the preparation of a nanoparticle suspension or a vesicle dispersion in an aqueous medium consists of three steps: a) assembly of the Nanocarriers (for example by injections, film hydration, or reverse phase evaporation), b) purification (exemplary: chromatography, dialysis or ultrafiltration), and c) concentration like ultrafiltration or evaporation.

This short review provides examples of recent literature dealing with the preparation of Nanocarriers. Particular focus is laid on the concentration and purification steps which were performed via ultrafiltration with Sartorius Vivaspin[®] or Vivaflow[®] devices with different pore sizes (respectively molecular weight cut-off, MWCO). The Vivaspin[®] portfolio spans a volume range from 0.5 mL to up to 20 mL, whereas the Vivaflow[®] system covers volumes from 0.05 liters to up to 5 liters. Thus, Sartorius offers an unrivalled wide range of processable sample volumes, membrane materials and MWCOs to meet the different requirements of their intended use. Challenges in this context are buffer exchange after synthesis, desalting and washing,^{10,11} exclusion of solubilized compounds,^{12,13,14} or aggregates.¹⁵

Purification is essential to obtain isosmotic conditions for in vivo applications to prevent aggregation or agglomeration and to remove free toxic drugs, ligands, or other substrates potentially triggering side effects. Concentration steps are essential to adjust the amount of pharmaceutical active ingredient in the drug to achieve the anticipated therapeutic or diagnostic effect.

During purification, the separation of free substances (starting material) from the desired Nanocarriers via size-exclusion chromatography (SEC) leads to an unavoidable dilution and to the necessity of a subsequent concentration step. In contrast, diafiltration purifies without significant dilution but a concentration step can still be mandatory, if higher Nanocarrier concentrations are necessary. Both separation methods require a quite extensive costly and time-consuming manual handling. This drawback is overcome by the ultrafiltration utilized by centrifugation in Vivaspin[®] or with a peristaltic pump for the Vivaflow[®] System. This technique is less expensive and quickly performed with very little manual input. Noteworthy is that purification and concentration steps are performed simultaneously.¹⁶

After the Nanocarrier is purified the determination of drug loading (conjugation or encapsulation efficiency) is commonly performed. The conjugation or encapsulation efficiency is one of the reference values to describe and characterize Nanocarriers. Other important properties are the zeta potential and the size distribution determined via photon correlation spectroscopy (PCS), high-resolution transmission electron microscopy (HRTEM) imaging, or via dynamic light scattering (DLS). Prior to performing these different characterizations a successful purification and concentration of the suspension or dispersion is essential.

In the following table you can find an overview of publications using ultrafiltration steps for the purification and concentration of different kinds of Nanocarriers. This table will also give you a guidance on which MWCOs to use.



* In German "Zauberkegel", opera "Freischütz" by Carl Maria von Weber

Table 1 summarizes examples of Nanocarrier ultrafiltration applications with Sartorius Vivaspin® or Vivaflow®:

Nanocarrier: Nanoparticle, Vesicle, Micelle	Size distribution obtained via (HR) TEM or DLS, Z-Average via PCS and others-if reported	Application	Sartorius Ultrafiltration Device	MWCO	Ultrafiltration purpose	Ref.
Nanoparticles from metal, metal oxides and functionalized metals						
Iron oxides nanoparticles with cisplatin-bearing polymer coating	SD: 4.5 ± 0.9 nm via X-Ray-Diffraction Analysis	Magnetic resonance imaging	Vivaspin® 20	100 kDa	Purification and concentration step	7
Functionalized iron oxide nanoparticles	SD: 38 and 40 nm via DLS	Stem cell gene therapy and tracking	Vivaspin® 20	100 kDa	Washing step	8
Gold nanoparticles	SD: 0.8 – 10.4 nm via Atomic Force Microscopy	Antimicrobial activity	Vivaspin® 20	5 kDa	Purification step	17
Protein coated gold nanoparticles	SD: 15 and 80 nm via TEM	Drug delivery	Vivaspin® 6	10 kDa	Separation of Nanoparticles Dyes and Washing	18
Functionalized gold nanoparticles	Core-SD: 2 nm via TEM	Targeted imaging tool and antigen delivery	Vivaspin®	10 kDa	Purification step	19
Functionalized gadolinium-based nanoparticles	Z-Average: 1.1 ± 0.6 nm and 4 – 14 nm	Diagnostic and therapeutic application	Vivaspin®	5 kDa, 10 kDa	Purification and Concentration	20, 21
Functionalized nanocrystals	SD: 10 to 20 nm	Quantum dots for imaging	Vivaspin®	300 kDa and 50 kDa	Separation of quantum dots-antibody conjugates from starting material (prior to enumeration)	9
Nanoparticles from polymers, functionalized polymers and polymersomes						
Polymer based Nanoparticles		Drug delivery	Vivaspin®	30 kDa	Purification and Concentration	22
Curdlan coated polymer nanoparticles	Z-Average: 280 – 480 nm depending on the composition	Macrophage stimulant activity and drug delivery	Vivaspin® 20	3 kDa	Washing	23
Docetaxel-carboxymethylcellulose Polymer Nanoparticles	Z-Average: 118 ± 1.8 nm	Anti-cancer efficacy studies	Vivaspin®	10 kDa	Concentration step	24
Functionalized Polymersomes	Z-Average: 185 nm	Surface functionalization studies	Vivaspin® 20	10 kDa	Concentration step	3
Lipid Nanoparticles and Liposomes						
Liposomes and micelles	Z-Average: 100 nm for Liposomes and 15 nm for micelles	Ischemia-reperfusion injury	Vivaspin® 20	100 kDa	Concentration step	25
Extracellular vesicles (Exosomes and microvesicles)	Exosomes: 70 – 150 nm Microvesicles: 100 – 1000 nm	Paper provides a general protocol	Vivaflow® 50R	100 kDa	Diafiltration and Concentration	26
Bacterial outer membrane vesicles	SD: 124 nm via TRPS	Tunable resistive pulse sensing (TRPS) Analysis	Vivaflow® 200	100 kDa	Buffer exchange and concentration step	27
Bacterial outer membrane vesicles		Basic research	Vivaspin® 20 and 500	100 kDa	Buffer exchange and concentration step	28
Bacterial outer membrane vesicles	SD: 95 nm	Basic research	Vivaflow® 200	100 kDa	Buffer exchange and concentration step	29
Bacterial outer membrane vesicles	SD: 50 – 150 nm via TEM	Basic research	Vivaspin®	100 kDa	Buffer exchange and concentration step	30
Liposomes		Drug delivery	Vivaspin®	100 kDa	External buffer exchange	2
Urinary exosomes	size of exosomes <100 nm	Preparation of urinary exosomes	Vivaspin® 20 and 500	100 kDa	Concentration	31
Micelles						
Micelles		Drug delivery	Vivaspin®	30 kDa	Separation of free substrate and concentration step	4
Hydrophobic drug micelles based on polymers	SD via DLS: 39 – 165 nm depending on compound in use	Drug delivery	Vivaflow®		Removal surfactant	14
Protein Nanoparticles						
Protein Nanoparticles	SD: 20 – 40 nm via DLS	Drug carrier studies	Vivaspin® 500	3 kDa	Separation of the free from the encapsulated drug (Drug binding quantification by subsequent UV Vis analysis)	32

SD = Size distribution

References

1. Strebhardt, K. & Ullrich, A.: Paul Ehrlich's magic bullet concept: 100 years of progress. *8*, 473–480 (2008).
2. Jakoby, J., Beuschlein, F., Mentz, S., Hantel, C. & Süss, R.: Liposomal doxorubicin for active targeting: Surface modification of the nanocarrier evaluated in vitro and in vivo – challenges and prospects. *Oncotarget* *6*, (2015).
3. Klermund, L., Poschenrieder, S. T. & Castiglione, K.: Simple surface functionalization of polymersomes using non-antibacterial peptide anchors. *J. Nanobiotechnology* *14*, 48 (2016).
4. Mulder, W. J. M. et al.: Molecular imaging of macrophages in atherosclerotic plaques using bimodal PEG-micelles. *Magn. Reson. Med.* *58*, 1164–1170 (2007).
5. Murthy, S. K.: Nanoparticles in modern medicine: state of the art and future challenges. *Int. J. Nanomedicine* *2*, 129–41 (2007).
6. Voigt, R. & Fahr, A.: *Pharmazeutische Technologie für Studium und Beruf.* Deutscher Apotheker Verlag, 10th Edition (2010).
7. Unterweger, H. et al.: Development and characterization of magnetic iron oxide nanoparticles with a cisplatin-bearing polymer coating for targeted drug delivery. *Int J Nanomedicine* *9*, 3659–3676 (2014).
8. Park, W. et al.: Multi-modal transfection agent based on monodisperse magnetic nanoparticles for stem cell gene delivery and tracking. *Biomaterials* *35*, 7239–7247 (2014).
9. Chalmers, N. I. et al.: Use of quantum dot luminescent probes to achieve single-cell resolution of human oral bacteria in biofilms. *Appl. Environ. Microbiol.* *73*, 630–636 (2007).
10. Hoffman, L. W., Andersson, G. G., Sharma, A., Clarke, S. R. & Voelcker, N. H.: New insights into the structure of PAMAM dendrimer/gold nanoparticle nanocomposites. *Langmuir* *27*, 6759–6767 (2011).
11. Rademacher, T. & Williams, P.: Nanoparticle-peptide compositions. (2014).
12. Allard, E. & Larpent, C.: Core-shell type dually fluorescent polymer nanoparticles for ratiometric pH-sensing. *J. Polym. Sci. Part A Polym. Chem.* *46*, 6206–6213 (2008).
13. Prach, M., Stone, V. & Proudfoot, L.: Zinc oxide nanoparticles and monocytes: Impact of size, charge and solubility on activation status. *Toxicol. Appl. Pharmacol.* *266*, 19–26 (2013).
14. Zhang, Y. et al.: Therapeutic surfactant-stripped frozen micelles. *Nat Commun* *7*, 11649 (2016).
15. Klasson, A. et al.: Positive MRI contrast enhancement in THP-1 cells with Gd2O3 nanoparticles. *Contrast Media Mol. Imaging* *3*, 106–111 (2008).
16. Simonoska Crcarevska, M. et al.: Definition of formulation design space, in vitro bioactivity and in vivo biodistribution for hydrophilic drug loaded PLGA/PEO-PPO-PEO nanoparticles using OFAT experiments. *Eur. J. Pharm. Sci.* *49*, 65–80 (2013).
17. Boda, S. K. et al.: Cytotoxicity of Ultrasmall Gold Nanoparticles on Planktonic and Biofilm Encapsulated Gram-Positive Staphylococci. *Small* *11*, 3183–3193 (2015).
18. Schäffler, M. et al.: Blood protein coating of gold nanoparticles as potential tool for organ targeting. *Biomaterials* *35*, 3455–3466 (2014).
19. Arosio, D. et al.: Effective targeting of DC-sign by α -fucosylamide functionalized gold nanoparticles. *Bioconjug. Chem.* *25*, 2244–2251 (2014).
20. Miladi, I. et al.: Biodistribution of ultra small gadolinium-based nanoparticles as theranostic agent: application to brain tumors. *J. Biomater. Appl.* *28*, 385–94 (2013).
21. Faure, A. C. et al.: Control of the in vivo biodistribution of hybrid nanoparticles with different poly(ethylene glycol) coatings. *Small* *5*, 2565–2575 (2009).
22. Benita, S., Debotton, N. & Goldstein, D.: Nanoparticles for Targeted Delivery of Active Agent. (2008).
23. Tukulula, M. et al.: Curdlan-conjugated PLGA nanoparticles possess macrophage stimulant activity and drug delivery capabilities. *Pharm. Res.* *32*, 2713–2726 (2015).
24. Ernsting, M. J., Tang, W. L., MacCallum, N. W. & Li, S. D.: Preclinical pharmacokinetic, biodistribution, and anti-cancer efficacy studies of a docetaxel-carboxymethylcellulose nanoparticle in mouse models. *Biomaterials* *33*, 1445–1454 (2012).
25. Geelen, T., Paulis, L. E., Coolen, B. F., Nicolay, K. & Strijkers, G. J.: Passive targeting of lipid-based nanoparticles to mouse cardiac ischemia-reperfusion injury. *Contrast Media Mol. Imaging* *8*, 117–126 (2013).
26. Corso, G. et al. Reproducible and scalable purification of extracellular vesicles using combined bind-elute and size exclusion chromatography *Scientific Rep.* *7*, 11561 (2017).
27. Bogomolny, E. et al.: Analysis of bacteria-derived outer membrane vesicles using tunable resistive pulse sensing. *Prog. Biomed. Opt. Imaging – Proc. SPIE* *9338*, 4–9 (2015).
28. Blenkiron, C. et al.: Uropathogenic Escherichia coli releases extracellular vesicles that are associated with RNA. *PLoS One* *11*, 1–16 (2016).
29. Twu, O. et al.: Trichomonas vaginalis Exosomes Deliver Cargo to Host Cells and Mediate Host:Parasite Interactions. *PLoS Pathog.* *9*, 22–24 (2013).
30. Tong, T. T., Mörgelin, M., Forsgren, A. & Riesbeck, K.: Haemophilus influenzae Survival during Complement-Mediated Attacks Is Promoted by Moraxella catarrhalis Outer Membrane Vesicles. *J. Infect. Dis.* *195*, 1661–1670 (2007).
31. Cheruvanky, A. et al. Rapid isolation of urinary exosomal biomarkers using a nanomembrane ultrafiltration concentrator. *Am. J. Physiol. Renal. Physiol.* *292*, F1657–F1661, (2007).
32. Achilli, E. et al.: Preparation of protein nanoparticle by dynamic aggregation and ionizing-induced crosslinking. *Colloids Surfaces A Physicochem. Eng. Asp.* *486*, 161–171 (2015).

Sartorius Lab Instruments
GmbH & Co. KG
Otto-Brenner-Strasse 20
37079 Goettingen, Germany
Phone +49.551.308.0
email: LF-info@sartorius.com
www.sartorius.com

USA Toll-free +1.800.635.2906
UK +44.1372.737159
France +33.1.70.62.50.00
Italy +39.0362.5557.11
Spain +34.913.586.095
Russian Federation +7.812.327.53.27
Japan +81.3.3740.5408

Specifications subject to change without notice.
Copyright Sartorius Lab Instruments GmbH & Co. KG.
Printed in the EU on paper bleached without chlorine.
Publication No.: SL-4001-e180502
Order No.: 85037-559-29
Ver. 05 | 2018