



Evaluation of a Centrifugal Filtration Device in the Measurement of Free Testosterone

July 29, 2015

Elizabeth Olson, BS MT(ASCP)

Sihe Wang, PhD DABCC FACB

Rob Kreller, BS MT(ASCP)


Michelle Strizzi, BS

Cleveland Clinic Background

- Testosterone
 - 2014 annual tests - 28,200
 - Automation (Enzyme Immunoassay)
- Free testosterone
 - 2014 annual tests - 21,000
 - Indirect measurement
 - Manual (Radiolimmunoassay)
 - Ultrafiltration
 - Millipore Centrifree™ Ultrafiltration device




Testosterone

	Male		Female
Produced	<ul style="list-style-type: none">• Testes• Adrenal glands		<ul style="list-style-type: none">• Ovaries• Adrenal glands
Responsibility	<ul style="list-style-type: none">• Development of external genitalia• Secondary sex characteristics		<ul style="list-style-type: none">• Estrogen precursor• Hormonal balance
Why	<ul style="list-style-type: none">• Infertility (low sperm counts)• Sexual problems (decreased sex drive/erectile dysfunction)• Premature/delayed puberty• Osteoporosis• Monitoring of testosterone therapy		<ul style="list-style-type: none">• Irregular menstrual cycles• Development of male features (hirsutism)

Total vs Free Testosterone

- Total testosterone = bound + free testosterone
- Majority of testosterone is bound
 - Sex hormone binding globulin (SHBG)
 - Albumin
- Free Testosterone
 - Unbound



Free Testosterone

- Metabolically active
- SHBG driven
 - Increase in SHBG = Decrease in Free Testosterone
 - Decrease in SHBG = Increase in Free Testosterone
- Testing importance
 - Androgen deficiency in males
 - Hypogonadism
 - Androgen excess in females
 - Polycystic Ovary Syndrome
 - Hirsutism

Free Testosterone Separation

Two methods of separation

- Equilibrium Dialysis (gold standard)
 - Time consuming
 - Technically difficult
 - Error prone
 - Not suitable for high volume testing
- Ultrafiltration
 - Faster
 - Technically easier
 - Centrifuge limitation




Ultrafiltration Device Evaluation

- Centrifree™ (EMD Millipore)
 - Current ultrafiltration device in use
 - No quality/testing issues with the device
- Vivafree™ (VivaProducts)
 - Proposed ultrafiltration device
 - Filter based on molecular weight of 30,000 daltons



Mechanism




Free Testosterone Comparison


- 53 serum samples
- Filtrates measured by Radiolimmunoassay

Sample Size	300 µL
Centrifuge Speed	4500 rpm
Centrifuge Time	15 min
Filtrate sample size	100 µL

Comparison Results -15 min spin



Slope = 0.963
Corr Coef (R) = 0.9685



Average % Bias = -9.71

Intraday Precision -15 min spin


Bio-Rad controls assayed

- High control
 - Samples tested (n) = 4
 - Mean 141.75 pg/mL
 - Observed CV = 5.8%
- Low Control
 - Samples tested (n) = 4
 - Mean 12.08 pg/mL
 - Observed CV = 3.0%


Conclusion

- High negative bias (-9.71%)
- Repeat comparison
 - Adjust spin time of Vivafree™ device to 20 min
 - 39 serum samples

Comparison Results - 20 min spin



Slope = 1.047
Corr Coef (R) = 0.9748



Average % Bias = 5.78

Intraday Precision -20 min spin

Bio-Rad controls assayed

- High control
 - Samples tested (n) = 10
 - Mean 188.09 pg/mL
 - Observed CV = 6.5%
- Low Control
 - Samples tested (n) = 10
 - Mean 16.02 pg/mL
 - Observed CV = 8.3%

Total Precision – 20 min spin

Bio-Rad controls assayed

- High control
 - Samples tested (n) = 52
 - Mean = 144.3 pg/mL
 - CV = 12.56%
- Low control
 - Samples tested (n) = 53
 - Mean = 13.0 pg/mL
 - CV = 10.18%

Conclusion

- The VivaFree™ device from Viva Products is comparable to the Millipore device for the determination of Free Testosterone
 - 20 min centrifugation time
 - Slope 1.047
 - r 0.9748
 - Precision \leq 12.56%
- Substantial cost savings
 - 35% reduction
 - \$42,000 annually
 - Better value
 - No additional tech time

