



Proline® Plus Mechanical Pipettes

Dependable Durability



Fixed volume Proline® Plus pipette



Proline® Plus has a comfortable handle and an ergonomic finger support for effortless pipetting.



Safe-Cone Filters protect the pipette from contamination, and should be changed regularly.

The Proline® Plus mechanical pipette family is designed to offer comfort and quality for your everyday manual pipetting. It combines durable construction with ease and lightness of use, and is therefore the perfect choice for liquid-handling professionals and students alike. In addition, it has the widest pipette range, including fixed volume pipettes for when volumes need to be ready-set to avoid errors.

Ergonomic Design

Proline® Plus has low pipetting forces, a comfortable handle and an ergonomic finger support for effortless pipetting. Good fit in hand minimizes the grip force needed to hold the pipette, thereby reducing the risk of strain injury.

Loading Tips on Multi-channel Pipettes with Minimum Force

Loading and ejecting tips with multi-channel pipettes requires relatively high forces. Proline® Plus multi-channel pipettes have spring loaded tip cones – the Optiload mechanism. Optiload allows tips to be loaded and ejected with minimum force, which reduces the risk of hand injuries. Moreover, it secures even tip sealing onto every individual tip cone.

Protect the Pipette with Safe-Cone Filters

The replaceable Safe-Cone Filter located inside the tip-cone prevents aerosols and fluids from penetrating the pipette, also in case of over-aspiration. The use of Safe-Cone Filters lengthens the maintenance interval of the pipette. Safe-Cone Filters are available for all Proline® Plus models greater than 10 µL. They must be replaced regularly, and always in case of overaspiration.

Easy Maintenance and Adjustment

No opening tools are needed for cleaning and maintaining Proline® Plus pipettes, and only three parts need to be cleaned. These pipettes are also easy to adjust using the adjustment tool that is provided with the pipette.

Features and Benefits



Proline® Plus pipettes have only three parts that need cleaning and maintenance.

- Low pipetting forces that prevent repetitive injury and improve results in long pipetting series
- Ergonomic finger support minimizes the grip force needed to hold the pipette
- Optiload mechanism in multi-channel models for easy and light tip loading with perfect tip sealing as well as 5,000 and 10,000 μL single channel pipettes
- Wide variety of adjustable single and multi-channel models as well as fixed single-channel models
- Volume range from 3 μL (5 μL for fixed) to 10 mL
- Easy volume setting with click stop mechanism
- Big volume display
- Color-coding of different volumes to ease the selection of matching pipette tips
- Safe-Cone Filters available for models > 10 μL
- Fully autoclavable without disassembly
- Simple to clean and maintain with only three parts to disassemble
- Materials have high chemical and UV-resistance to secure long life-time of the pipette



Ordering Information

Proline® Plus Variable Volume Mechanical Pipettes

Tip Selection Guide

Proline® Plus Order Code	Channels	Volume Range	Increment	Pipettes and Tips	Optifit Tip	Safetyspace® Tip	Safe-Cone Filters	
		(µL)	(µL)	Color-Code	Volume	Volume	Standard	Plus
728010	1	0.3 – 3	0.002	■	10 µL	10 µL	–	–
728020	1	1 – 10	0.01	■	10 µL	10 µL	–	–
728030	1	2 – 20	0.02	■	200 µL	20 µL	721014	–
728040	1	5 – 50	0.10	■	200 µL	120 µL	721008	721018
728050	1	10 – 100	0.10	■	200 µL	120 µL	721008	721018
728060	1	20 – 200	0.20	■	200 µL	200 µL	721007	721017
728070	1	100 – 1,000	1.00	■	1,000 µL	1,000 µL	721006	721016
728080	1	500 – 5,000	10.0	■	5,000 µL	5,000 µL	721005	721015
728090	1	1,000 – 10,000	20.0	■	10,000 µL	–	721005	721015
728120	8	1 – 10	0.01	■	10 µL	10 µL	–	–
728220	12							
728130	8	10 – 100	0.10	■	200 µL	120 µL	721008	721018
728230	12							
728140	8	30 – 300	0.20	■	350 µL	300 µL	721007	721017
728240	12							

For pipette specifications (incl. Systematic and Random error limits), please see product datasheet, available at sartorius.com



Ordering Information

Tip Selection Guide

Proline® Plus FIXED Volume, Single Channel

Order Code	Channels	Volume	Pipettes and Tips	Optifit Tip	Safetyspace® Tip	Safe-Cone Filters	
		(µL)	Color-Code	Volume	Volume	Standard	Plus
728515	1	5	■	10 µL	10 µL	-	-
728520	1	10	■	10 µL	10 µL	-	-
728530	1	20	■	200 µL	20 µL	721014	-
728535	1	25	■	200 µL	120 µL	721008	721018
728545	1	50	■	200 µL	120 µL	721008	721018
728550	1	100	■	200 µL	120 µL	721008	721018
728560	1	200	■	200 µL	200 µL	721007	721017
728565	1	250	■	1,000 µL	1,000 µL	721006	721016
728567	1	500	■	1,000 µL	1,000 µL	721006	721016
728570	1	1,000	■	1,000 µL	1,000 µL	721006	721016
728575	1	2,000	■	5,000 µL	5,000 µL	721005	721015
728580	1	5,000	■	5,000 µL	5,000 µL	721005	721015
728590	1	10,000	■	10,000 µL		721005	721015

For pipette specifications (incl. Systematic and Random error limits), please see product datasheet, available at [sartorius.com](https://www.sartorius.com)



Proline® Plus Fixed Volume Pipettes



Proline® Plus Multipacks – Complete Sets of Pipettes and Accessories

Proline® Plus Multipacks offer are sets of mechanical pipettes in an affordable package, including a Linear Stand and racks of matching tips.

These Multipacks allow the affordable renewal of existing pipettes, or to set up new workstations, with a set of highly ergonomic pipettes with high accuracy and precision.

Ordering Information

Proline® Plus Mechanical Pipette Multipacks

Multipack Order Code	Proline® Plus Pipettes	Sartorius Optifit Tips	Accessories
Proline® Plus Pipette 3-pack 10 LH-728670	1-channel ▪ 1 – 10 µL ▪ 10 – 100 µL ▪ 100 – 1,000 µL	Tip Racks (96 tips) ▪ 10 µL ▪ 200 µL ▪ 1,000 µL	Linear Stand
Proline® Plus Pipette 3-pack 20 LH-728671	1-channel ▪ 2 – 20 µL ▪ 20 – 200 µL ▪ 100 – 1,000 µL	Tip Racks (96 tips) ▪ 200 µL ▪ 1,000 µL	Linear Stand
Proline® Plus Pipette 4-pack 20 LH-728672	1-channel ▪ 1 – 10 µL ▪ 2 – 20 µL ▪ 20 – 200 µL ▪ 100 – 1,000 µL	Tip Racks (96 tips) ▪ 10 µL ▪ 200 µL ▪ 1,000 µL	Linear Stand
Proline® Plus Pipette 4-pack 100 LH-728673	1-channel ▪ 1 – 10 µL ▪ 10 – 100 µL ▪ 20 – 200 µL ▪ 100 – 1,000 µL	Tip Racks (96 tips) ▪ 10 µL ▪ 200 µL ▪ 1,000 µL	Linear Stand
Proline® Plus Pipette 5-pack 10 LH-728674	1-channel ▪ 1 – 10 µL ▪ 10 – 100 µL ▪ 20 – 200 µL ▪ 100 – 1,000 µL ▪ 500 – 5,000 µL	Tip Racks (96 tips) ▪ 10 µL ▪ 200 µL ▪ 1,000 µL Tip Rack (50 tips) ▪ 5,000 µL	Linear Stand
Proline® Plus Pipette 5-pack 20 LH-728675	1-channel ▪ 2 – 20 µL ▪ 10 – 100 µL ▪ 20 – 200 µL ▪ 100 – 1,000 µL ▪ 500 – 5,000 µL	Tip Racks (96 tips) ▪ 200 µL ▪ 1,000 µL Tip Rack (50 tips) ▪ 5,000 µL	Linear Stand
Proline® Plus Pipette 3+1-pack LH-728676	1-channel ▪ 1 – 10 µL ▪ 20 – 200 µL ▪ 100 – 1,000 µL 8-channel ▪ 30 – 300 µL	Tip Racks (96 tips) ▪ 10 µL ▪ 200 µL ▪ 350 µL ▪ 1,000 µL	Linear Stand
Proline® Plus Pipette 3-pack max LH-728677	1-channel ▪ 100 – 1,000 µL ▪ 500 – 5,000 µL ▪ 1,000 – 10,000 µL	Tip Racks ▪ 1,000 µL (96 tips) ▪ 5,000 µL (50 tips) ▪ 10,000 µL (35 tips)	Linear Stand